


Americot® 100% Cotton Neurosurgical Pattie



- High pliability and unmatched softness with only a single barium line
- Precision cut providing smooth edges, reducing the potential for loose fibers
- Made with 100% U.S.A. Cotton
- X-ray detectable with locator string

Unmatched Softness

For Delicate Tissue

Americot® neurosurgical patties are made with 100% U.S.A. cotton. Americot® offers excellent absorption, softness and protection with the slow-drying properties of cotton compared to other synthetic based neuro patties on the market.



Superior Absorbency

100% U.S.A. Cotton

Americot® is capable of holding greater than 11 times its dry weight in fluid when fully saturated.



Product Codes & Sizes		
Americot code	size-inches	size-mm
20-01S	¼ x ¼	6 x 6
20-04S	½ x ½	13 x 13
20-05S	½ x 1	13 x 25
20-06S	½ x 1½	13 x 38
20-07S	½ x 3	13 x 76
20-08S	¾ x ¾	19 x 19
20-09S	1 x 1	25 x 25
20-10S	1 x 3	25 x 76
20-12S	3 x 3	76 x 76
20-13S	¾ x 3	19 x 76
20-14S	3/8 x ¾	10 x 19
20-19S	1/8 x 1/8	3 x 3
20-22S	3/8 x 3/8	10 x 10
20-25S	3/8 x 1	10 x 25
20-26S	3/8 x 2	10 x 51
20-27S	¾ x 2	19 x 51
20-30S	¼ x 1 ½	6 x 38
20-32S	½ x 2	13 x 51

Americot strips without locator string	size-inches	size-mm
21-01	¼ x 6	6 x 152
21-02	½ x 6	13 x 152
21-04	1 x 6	25 x 152

Manufactured in Salem, Massachusetts, U.S.A.
 Packaged: 200 per box (20 pouches of 10 patties)
 *All Americot® patty strings are blue except 20-01S and 20-19S, which are gold and 21-01, 21-02, and 21-01, 21-02, and 21-04, which are white.

Also Available in the UK from Minim Healthcare



Cotton Balls

X-Ray Detectable



Delicot® Neuro Pattie

Ultra Thin for Microneurosurgery



Ray-Cot® Neuro Pattie

Absorbent and Pliable Cottonoid