


X-Ray Detectable Strung Cotton Balls



- Each sterile strung cotton ball contains an x-ray detectable marker for increased safety and accountability.
- Available with or without locator strings. The locator strings allow for easier traceability post-procedure.
- Available in 8 sizes ranging from 3mm to 51mm, this eliminates the need to spend time manually tearing tufts of cotton during procedures which can impact accountability and lead to gossypiboma or textiloma

Double Strung

For Increased Safety

Cotton balls are x-ray detectable reducing the risk of gossypiboma or longer operating times. Commonly used in a variety of surgical procedures where a pattie can't be used.



Product Codes & Sizes

Strung cotton balls code	size-inches	size-mm
30-07 / x-small	1/8	3
30-00 / small	1/4	6
30-03 / medium	3/4	19
30-04 / large	1	25
30-05 x-large	1 1/4	32
Unstrung cotton balls code	size-inches	size-mm
31-07 / x-small	1/8	3
31-00 / small	1/4	6

Manufactured in Salem, Massachusetts, U.S.A.

100 per box (20 pouches of 5)

Also Available in the UK from Minim Healthcare




Americot® Neuro Pattie

Absorbent and Pliable Cottonoid



Delicot® Neuro Pattie

Ultra Thin for Microneurosurgery



Ray-Cot® Neuro Pattie

Absorbent and Pliable Cottonoid