


## Ray-Cot® Neurosurgical Pattie



- 12% greater absorbency and 2 times stronger than the leading standard patty competitor
- Single barium line allows for superior pliability and softness
- Smooth precise edges reducing the potential for loose fibers
- Thinner locator string for improved visibility
- X-ray detectable with locator string

## Improved Standard Pattie

Soft and Absorbent

Ray-Cot® is an improved version of the standard rayon neurosurgical patty with superior strength and performance.


Ray-Cot® neurosurgical patties utilize a special rayon material that has superior pliability and absorption, while providing softness and protection to delicate tissue. The Ray-Cot® material is shown to be up to 2x stronger than the leading competitor and demonstrates up to 12% greater absorption by weight per absorption testing.



## Clean Precise Edges

Reducing Loose Fibers

Clean precision cuts provide smooth edges and reduce loose fibers. The single barium line allows for superior flexibility and pliability. Ray-Cot® neuro patties are a staple in surgical suites.



Product Codes & Sizes			Product Codes & Sizes		
RayCot code	size-inches	size-mm	RayCot strips without locator string	size-inches	size-mm
60-02	¼ x ¼	6 x 6	61-01	¼ x 6	6 x 152
60-04	¼ x 3	6 x 76	61-02	½ x 6	13 x 152
60-05	3/8 x 3	10 x 76	61-03	¾ x 6	19 x 152
60-06	½ x ½	13 x 13	61-04	1 x 6	25 x 152
60-07	½ x 1	13 x 25	61-05	2 x 6	51 x 152
60-08	½ x 1 ½	13 x 38	61-06	3 x 6	76 x 152
60-09	½ x 3	13 x 76	61-07	3 ½ x 6	89 x 152
60-10	¾ x ¾	19 x 19	61-08	1 ½ x 6	38 x 152
60-12	1 x 1	25 x 25			
60-13	1 x 3	25 x 76			
60-14	1 ½ x 1 ½	38 x 38			
60-15	3 x 3	76 x 76			
60-16	¼ x 1 ½	6 x 38			
60-22	1 ½ x 3	38 x 76			
60-29	½ x 2	13 x 51			

Manufactured in Salem, Massachusetts, U.S.A.  
Packaged: 200 per box (20 pouches of 10 patties)

## Also Available in the UK from Minim Healthcare



Cotton Balls

X-Ray Detectable



Delicot® Neuro Pattie

Ultra Thin for Microneurosurgery



Americot® Neuro Pattie

Absorbent and Pliable Cottonoid