

Integra™ Omni-Tract®

Table Mounted
Upper Abdominal Retractor System

Limit uncertainty with
unparalleled upper abdominal
exposure with superior
subcostal margin elevation



INTEGRA™
LIMIT UNCERTAINTY



Integra™ Omni-Tract®

Table Mounted Upper Abdominal Retractor System

Table Mounted Retraction for Unparalleled Exposure in General Surgery Procedures


The Integra™ Omni-Tract® Upper Abdominal Retractor System, excellent retraction for obese patients and transplant cases

- Dual post frame for increased strength
- Angled side bars for low profile retraction
- Variety of retractor blades including sternal retractors
- Frame adjusts to fit operating room table extenders
- Fast setup and fast mid-surgery changes



Obesity Retractor Blade Kit:

- Optional obesity blade kit, deeper blades for large or obese patients for increased visibility
- Obesity blades are compatible with all Wishbone® style retractor systems




System can be used for Upper Abdominal surgeries including:

- Liver Transplants
- Repair or removal of the spleen, pancreas, or gallbladder
- Gastric Bypass
- Gastric resection/repair




Upper Abdominal Retractor System/Transplant System

Catalog Number	Qty	Recommended Set List
A. 4944	2	Omni-Clamp Sterile Field Post, Double Coupler, 24.3"
B. 4943	1	Standard Cross Bar
C. 4940	1	Hinged Cross Bar
D. 4040	7	SnapClamp Blade Holder
E. 3759	2	Angled Side Arm, 8" x 18"
F. 3091	1	Doherty-Harrington Retractor Blade
G. 3094	2	Sternal Retractor, 3" x 3"
H. 3536	2	Mayo Swivel Retractor, 3-1/2" x 2-3/4"
J. 3623	1	Blade Holder with 3626 Knob, 9"
K. 3872	1	Malleable Swivel Retractor, 2" x 6-1/2"
K. 3876	1	Malleable Swivel Retractor, 3" x 6-1/2"



Obesity Blade Kit

Catalog Number	Qty	Recommended Set List
A. 3529	2	Splanchnic Retractor, 2" x 8"
B. 4929	2	Renal Vein Swivel Retractor, 1" x 8"
C. 3531	1	Deaver Retractor, 2-3/4" x 8"
D. 3536	2	Mayo Swivel Retractor, 3-1/2" x 2-3/4"
E. 3856	1	Fence Swivel Retractor, 4" x 8"
F. 3874	1	Malleable Swivel Retractor, 2" x 8"
G. 3878	1	Malleable Swivel Retractor, 3" x 8"
H. 3882	1	Malleable Swivel Retractor, 4" x 8"



Integra™ Omni-Tract®

Table Mounted Upper Abdominal Retractor System

Recommended Containers

Catalog Number	Description
3424	Pan/Cover Set, 12-1/4"x18-1/2"x3-5/8"
3422	Pan/Cover Set, 12-1/2"x26-1/2"x3-5/8"
730-312	Full Size Container, Perforated Bottom, 23" x 11" x 6"
730-302	Full Size Container, Solid Bottom, 23" x 11" x 6"
750-302	Full Size Wire/Perforated Basket, 4"

*Must order to accommodate the 4944 tall posts



Pan/Cover
*3422



Pan/Cover
3424



Full Size Container, Solid Bottom
730-312 and 730-302



Full Size Wire/Perforated Basket
750-302

For more information or to place an order, please contact:
Integra ■ 4849 White Bear Parkway, St. Paul, MN 55110
800-367-8657 USA ■ 651-287-4300 outside USA ■ 651-287-4400 fax
integralife.com

