

Integra® Omni-Tract®

Table Mounted
Wishbone® Retractor Systems

Limit uncertainty with
flexibility for any approach



INTEGRA®
LIMIT UNCERTAINTY

Integra® Omni-Tract®



Table Mounted Wishbone® Retractor Systems

Limit uncertainty with flexibility for any approach.

Manufacturer of the Omni-Tract® Wishbone® Table Mounted Retractor Systems

Omni-Tract Surgical has over 30 years of experience in table mounted retractor design. The concept of securing the retractor frame to the O.R. table rail allows you to place retractor blades at any angle with exact pressure for unilateral, bilateral, single or multi-level exposures.

Table mounted retraction eliminates hand held retractors and self retaining retractor migration. Wishbone® retractor systems perform consistently on all patient sizes, and can be used in multiple surgical specialties.




Retract from every direction by simply adding our hinged extension arm.

Table Mounted Wishbone® Retractor Systems


General Surgery	3-4
Vascular Surgery	5-6
Urologic Surgery	7-8
Pediatric Surgery	9-10
Containers	11

Table Mounted Retraction for Unparalleled Exposure in General Surgery Procedures



One-piece Wishbone® Retractor Systems

- One piece Wishbone® Frame allows body wall and deep retraction regardless of patient size and anatomical structure
- Multi positional blades firmly retract from any plane
- Low profile system follows contour of operative field
- Optional hinged extension arm encloses the operative site for 360° retraction
- Wishbone® Frames have SnapClamp™ clamps for quick blade placement.




Flexible Wishbone® Retractor Systems

- Flexible Wishbone® Frames have interchangeable, removable frame arms that allow you to maximize your surgical field
- Strength – no matter your patient size the Flexible frame arms are designed to remain stable, so you will always get the exposure you need
- Hinged extension arm is also available so you retract from every direction if needed

Wishbone® Systems can be used in general surgeries including:
Gastric Resection, Exploratory Laparotomy, Open Cholecystectomy, Splenectomy and Gastic Bypass procedures.

One Piece Wishbone® General Surgery System

Catalog Number	Qty	Description
A. 4080	1	Omni Clamp® Sterile Field Post
B. 4020	1	Standard Wishbone® Frame
C. 4040	8	SnapClamp™ Blade Holder
D. 3588	1	Hinged Extension Arm, 19"
E. 3335	2	Mayo Swivel Retractor, 2 3/4" x 2"
F. 3772	1	Fence Swivel Retractor, 4" x 6"
G. 3774	1	Splanchnic Swivel Retractor, 2" x 6"
H. 3784	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2"
I. 3788	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2"
J. 3790	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2"
K. 3794	1	Harrington Sweetheart Swivel Retractor, 2" x 5"
L. 3862	1	Kelly Modified Swivel Retractor, 3" x 2 1/2"
M. 3864	1	Kelly Modified Swivel Retractor, 3" x 4"



Flexible Wishbone® General Surgery System

Catalog Number	Qty	Description
A. 10244	1	Sterile Field Post
B. 11050	1	Support Arm
C. 10269	2	Standard Wishbone® Frame Arm
D. 10270	2	Small Wishbone® Frame Arm
E. 4040	8	SnapClamp™ Blade Holder
F. 11121	1	Hinged Extension Arm, 19"
G. 3335	2	Mayo Swivel Retractor, 2 3/4" x 2"
H. 3772	1	Fence Swivel Retractor, 4" x 6"
I. 3774	1	Splanchnic Swivel Retractor, 2" x 6"
J. 3784	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2"
K. 3788	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2"
L. 3790	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2"
M. 3794	1	Harrington Sweetheart Swivel Retractor, 2" x 5"
N. 3862	1	Kelly Modified Swivel Retractor, 3" x 2 1/2"
O. 3864	1	Kelly Modified Swivel Retractor, 3" x 4"






Table Mounted Retraction for Multi plane, Multi positional Hands Free Exposure in Vascular Surgery Procedure




One piece Wishbone® Retractor Systems

- One piece Wishbone® Frame allows body wall and deep retraction regardless of patient size and anatomical structure
- Multi positional blades firmly retract from any plane
- Low profile system follows contour of operative field
- Optional hinged extension arm encloses the operative site for 360° retraction
- Wishbone® Frames have SnapClamp™ clamps for quick blade placement.



Flexible Wishbone® Retractor Systems


- Flexible Wishbone® Frames have interchangeable, removable frame arms that allow you to maximize your surgical field
- Strength – no matter your patient size the Flexible frame arms are designed to remain stable, so you will always get the exposure you need
- Hinged extension arm is also available so you retract from every direction if needed



Wishbone® Systems can be used in vascular surgeries including: Abdominal Aortic Aneurysm and Iliofemoral Bypass procedures, and for small site procedures such as Carotid Endarterectomy, Femoral Popliteal and Thoracic Outlet.


One piece Wishbone® Vascular Surgery System

Catalog Number	Qty	Description
A. 4080	1	Omni Clamp® Sterile Field Post
B. 4020	1	Standard Wishbone® Frame
C. 4040	8	SnapClamp™ Blade Holder
D. 3588	1	Hinged Extension Arm, 19"
E. 3335	2	Mayo Swivel Retractor, 2 3/4" x 2"
F. 3397	1	Mayo Swivel Retractor, 3 1/2" x 2"
G. 3772	1	Fence Swivel Retractor, 4" x 6"
H. 3888	1	Splanchnic Swivel Retractor, 2" x 4 1/2"
I. 3774	1	Splanchnic Swivel Retractor, 2" x 6"
J. 4928	2	Renal Vein Swivel Retractor, 1" x 6"
K. 3784	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2"
L. 3788	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2"
M. 3790	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2"
N. 3794	1	Harrington Sweetheart Swivel Retractor, 2" x 5"




Flexible Wishbone® Vascular Surgery System

Catalog Number	Qty	Description
A. 10244	1	Sterile Field Post
B. 11050	1	Support Arm
C. 10269	2	Standard Wishbone® Frame Arm
D. 10270	2	Small Wishbone® Frame Arm
E. 4040	8	SnapClamp™ Blade Holder
F. 11121	1	Hinged Extension Arm, 19"
G. 3335	2	Mayo Swivel Retractor, 2 3/4" x 2"
H. 3397	1	Mayo Swivel Retractor, 3 1/2" x 2"
I. 3772	1	Fence Swivel Retractor, 4" x 6"
J. 3888	1	Splanchnic Swivel Retractor, 2" x 4 1/2"
K. 3774	1	Splanchnic Swivel Retractor, 2" x 6"
L. 4928	2	Renal Vein Swivel Retractor, 1" x 6"
M. 3784	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2"
N. 3788	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2"
O. 3790	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2"
P. 3794	1	Harrington Sweetheart Swivel Retractor, 2" x 5"




Flexible Wishbone® Stoney Peripheral Vascular Surgery System

Catalog Number	Qty	Description
A. 10244	1	Sterile Field Post
B. 11050	1	Support Arm
C. 10270	2	Small Wishbone® Frame Arm
D. 4040	4	SnapClamp™ Blade Holder
E. 4958	1	Rake Retractor, 7/8" x 1/2"
F. 4957	1	Rake Retractor, 7/8" x 1"
G. 4956	1	Richardson Stoney Ret. Micro-Adjust, 3/4" x 2"
H. 4959	1	Richardson Swivel Retractor, 3/4" x 1"
I. 4960	1	Richardson Swivel Retractor, 1" x 1"
J. 4961	1	Richardson Swivel Retractor, 1" x 2"




*Also available with the one piece Wishbone® Frame

Table Mounted Retraction for Unparalleled Exposure in Urologic Surgery Procedures




One-piece Wishbone® Retractor Systems

- One piece Wishbone® Frame allows body wall and deep retraction regardless of patient size and anatomical structure
- Multi positional blades firmly retract from any plane
- Optional hinged extension arm encloses the operative site for 360° retraction field
- Wishbone® Frames have SnapClamp™ clamps for quick blade placement.



Flexible Wishbone® Retractor System


- Flexible Wishbone® Frames have interchangeable, removable frame arms that allow you to maximize your surgical field
- Strength – no matter your patient size the Flexible frame arms are designed to remain stable, so you will always get the exposure you need
- Hinged extension arm is also available so you retract from every direction if needed



Wishbone® Systems can be used in urologic surgeries including: Radial Retropubic Prostatectomy, Nephrectomy and Cystectomy procedures, and GYN procedures such as Perineal Prostatectomy, Bladder Suspension and Vaginal Hysterectomy.


One piece Wishbone® Urologic Surgery System

Catalog Number	Qty	Description
A. 4080	1	Omni Clamp® Sterile Field Post
B. 4020	1	Standard Wishbone® Frame
C. 4040	8	SnapClamp™ Blade Holder
D. 3588	1	Hinged Extension Arm, 19"
E. 3883	1	Mayo Swivel Retractor, 2" x 2"
F. 3335	2	Mayo Swivel Retractor, 2 3/4" x 2"
G. 3772	1	Fence Swivel Retractor, 4" x 6"
H. 3784	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2"
I. 3854	1	Deaver Swivel Retractor, 1 1/2" x 6"
J. 3788	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2"
K. 3790	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2"
L. 3792	1	Marshall Bladder Malleable Swivel Retractor, 3" x 6"
M. 3794	1	Harrington Sweetheart Swivel Retractor, 2" x 5"
N. 3796	1	Harrington Sweetheart Swivel Retractor, 2 1/2" x 7"
O. 3798	1	Splanchnic Swivel Retractor, 1 1/2" x 4"




Flexible Wishbone® Urologic Surgery System

Catalog Number	Qty	Description
A. 10244	1	Sterile Field Post
B. 11050	1	Support Arm
C. 10269	2	Standard Wishbone® Frame Arm
D. 10270	2	Small Wishbone® Frame Arm
E. 4040	8	SnapClamp™ Blade Holder
F. 11121	1	Hinged Extension Arm, 19"
G. 3883	1	Mayo Swivel Retractor, 2" x 2"
H. 3335	2	Mayo Swivel Retractor, 2 3/4" x 2"
I. 3772	1	Fence Swivel Retractor, 4" x 6"
J. 3784	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2"
K. 3854	1	Deaver Swivel Retractor, 1 1/2" x 6"
L. 3788	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2"
M. 3790	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2"
N. 3792	1	Marshall Bladder Malleable Swivel Retractor, 3" x 6"
O. 3794	1	Harrington Sweetheart Swivel Retractor, 2" x 5"
P. 3796	1	Harrington Sweetheart Swivel Retractor, 2 1/2" x 7"
Q. 3798	1	Splanchnic Swivel Retractor, 1 1/2" x 4"




Flexible Wishbone® Pelvic Surgery System

Catalog Number	Qty	Description
A. 10244	1	Sterile Field Post
B. 11050	1	Support Arm
C. 10270	2	Small Wishbone® Frame Arm
D. 4040	5	SnapClamp™ Blade Holder
E. 11070	1	Extension Arm, 20"
F. 4977	1	Anterior Swivel Retractor, 1-1/2" x 3"
G. 4964	4	Malleable Swivel Retractor, 3/4" x 2"
H. 4965	4	Malleable Swivel Retractor, 3/4" x 3"
I. 4966	4	Malleable Swivel Retractor, 3/4" x 4"
J. 4975	1	Posterior Swivel Retractor, 1-3/4" x 3"
K. 4976	1	Posterior Swivel Retractor, 1-3/4" x 4"



*Also available with the one piece Wishbone® Frame


Table Mounted Pediatric Wishbone® Retractor Systems



One-piece Wishbone® Pediatric Retractor Systems

One piece Wishbone® Retractor system designed for pediatric procedures. Small, independently-adjustable, low-profile frame arms follow contours of the operative field

- Multi-plane, multi-positional hands-free retraction
- Accommodates variation in patient size
- Wishbone® Frames have SnapClamp™ clamps for quick blade placement.
- Can be used on infants, toddlers and adolescent patients



Flexible Wishbone® Pediatric Retractor System


Flexible Wishbone® Retractor System designed for pediatric patients of all sizes.

- Small, independently-adjustable and removable frame arms follow contours of any pediatric operative field
- Multi-plane, multi-positional hands-free retraction
- Easily accommodates the range of pediatric patients from infants to adolescents

Wishbone® Systems can be used in pediatric surgeries including:
Pediatric Transplant and Nissen Fundoplication procedures.


One Piece Wishbone® Pediatric Surgery System

Catalog Number	Qty	Description
A. 4080	1	Omni-Clamp® Sterile Field Post
B. 4985	1	Small Wishbone® Frame
C. 4040	6	SnapClamp™ Blade Holder
D. 4958	1	Rake Retractor, 7/8" x 1/2" x 6"
E. 4957	1	Rake Retractor, 7/8" x 1" x 6"
F. 4978	1	Mayo Swivel Retractor, 1" x 1" x 6"
G. 4979	2	Mayo Swivel Retractor, 2" x 1-1/2" x 6"
H. 4980	2	Mayo Swivel Retractor, 2" x 2" x 6"
I. 4981	2	Mayo Swivel Retractor, 2-3/4" x 2" x 6"
J. 4967	1	Malleable Swivel Retractor, 1" x 1" x 6"
K. 4963	1	Malleable Swivel Retractor, 1/2" x 1" x 6"
L. 4964	1	Malleable Swivel Retractor, 3/4" x 2" x 6"
M. 4965	1	Malleable Swivel Retractor, 3/4" x 3" x 6"
N. 4968	1	Malleable Swivel Retractor, 1" x 4" x 6"
O. 4969	1	Malleable Swivel Retractor, 1" x 5" x 6"
P. 4970	1	Malleable Swivel Retractor, 1" x 6" x 6"
Q. 4971	1	Malleable Swivel Retractor, 1-1/2" x 3" x 6"
R. 4972	1	Malleable Swivel Retractor, 1-1/2" x 4" x 6"
S. 4973	1	Malleable Swivel Retractor, 1-1/2" x 5" x 6"
T. 4974	1	Malleable Swivel Retractor, 1-1/2" x 6" x 6"
U. 4962	2	Splanchnic Swivel Retractor, 1-1/2" x 3-1/2" x 6"



Flexible Wishbone® Pediatric Surgery System

Catalog Number	Qty	Description
A. 10244	1	Sterile Field Post
B. 11050	1	Support Arm
C. 10270	2	Small Wishbone® Frame Arm
D. 4040	5	SnapClamp™ Blade Holder
E. 4958	1	Rake Retractor, 7/8" x 1/2" x 6"
F. 4957	1	Rake Retractor, 7/8" x 1" x 6"
G. 4978	1	Mayo Swivel Retractor, 1" x 1" x 6"
H. 4979	2	Mayo Swivel Retractor, 2" x 1-1/2" x 6"
I. 4980	2	Mayo Swivel Retractor, 2" x 2" x 6"
J. 4981	2	Mayo Swivel Retractor, 2-3/4" x 2" x 6"
K. 4967	1	Malleable Swivel Retractor, 1" x 1" x 6"
L. 4963	1	Malleable Swivel Retractor, 1/2" x 1" x 6"
M. 4964	1	Malleable Swivel Retractor, 3/4" x 2" x 6"
N. 4965	1	Malleable Swivel Retractor, 3/4" x 3" x 6"
O. 4968	1	Malleable Swivel Retractor, 1" x 4" x 6"
P. 4969	1	Malleable Swivel Retractor, 1" x 5" x 6"
Q. 4970	1	Malleable Swivel Retractor, 1" x 6" x 6"
R. 4971	1	Malleable Swivel Retractor, 1-1/2" x 3" x 6"
S. 4972	1	Malleable Swivel Retractor, 1-1/2" x 4" x 6"
T. 4973	1	Malleable Swivel Retractor, 1-1/2" x 5" x 6"
U. 4974	1	Malleable Swivel Retractor, 1-1/2" x 6" x 6"
V. 4962	2	Splanchnic Swivel Retractor, 1-1/2" x 3-1/2" x 6"



*Retractor blades and clamps are interchangeable with all Wishbone® frames

Integra® Omni-Tract®

Table Mounted Wishbone® Retractor Systems

Integra® Table Mounted Retractor Containers and Trays

Catalog Number	Description
----------------	-------------

Integra containers for TMR systems

730-612	Full Size Container, Perforated Bottom, 23" x 11" x 6"
730-602	Full Size Container, Solid Bottom, 23" x 11" x 6"
750-602	Full Size Wire/Perforated Basket, 4"
730-701	Wide Container, Solid Bottom, 120mm
730-711	Wide Container, Perforated Bottom, 120mm

TMR Trays

3424	Pan/Cover Set, 12¼" x 18½" x 3⅝"
3422R	Pan/Cover Set, 12½" x 26½" x 3⅝"



Pan/Cover

3424



Full Size Container, Solid Bottom

730-602



Full Size Wire/Perforated Basket

750-602

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-367-8657 ▪ 651-287-4400 fax

Outside USA 651-287-4300 ▪ 651-287-4400 fax

integralife.com

INTEGRA
LIMIT UNCERTAINTY