


Integra[®] MicroFrance[®]

Bouchayer Laryngoscopes

The Bouchayer Laryngoscopes offer a large proximal opening for binocular microscope vision, an elevated distal tip for retraction of the epiglottis and exposure of the anterior commissure, and an integrated light fiber for visualization.



Bouchayer Large Set

Reference	Description	Qty
MCL700	Bouchayer Laryngoscope Set, Large	1 Set
Set includes:		
MCL77	Bouchayer Laryngoscope, Large	1
MCL78	Light Fiber	1



Overall Length: 18cm
Tip Height (Retraction): 9cm
Proximal Opening: 42mm x 28mm

MCL77 Laryngoscope, Large

Bouchayer Medium Set

Reference	Description	Qty
MCL710	Bouchayer Laryngoscope Set, Medium	1 Set
Set includes:		
MCL79	Bouchayer Laryngoscope, Medium	1
MCL78	Light Fiber	1



Overall Length: 18cm
Tip Height (Retraction): 7.5cm
Proximal Opening: 39mm x 25mm

MCL79 Laryngoscope, Medium

Bouchayer Laryngoscope Accessories

Reference	Description	Qty
MCL76	Laryngoscope Suspension Apparatus Overall Length 35cm (13.8")	1
MCL35	Chest Plate for Suspension Apparatus 60mm Wide x 140mm Long	1
MCL78	Replacement Light Fiber for MCL77 and MCL79	1
3717023	Laryngoscope Instrument Tray 52 .5cm L x 25cm W x 17cm D Holds Two Laryngoscopes, Two Fiberoptic Cables, Suspension Apparatus, and Chest Plate	1



MCL76
Suspension Apparatus



3717023

The MicroFrance brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

Manufacturer:



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France