

Integra[®] MicroFrance[®]

Dessi Nasal Bipolar Forceps

Developed in conjunction with Patrick Dessi, MD

Hemostasis is one of the greatest challenges in endoscopic sinus surgery. The Dessi Nasal Bipolar Forceps feature a tip design that allows access to difficult areas of the nasal cavity, while the bipolar cautery is intended to remove tissue and control bleeding.



Integra® MicroFrance®

Dessi Nasal Bipolar Forceps

Dessi Bipolar Forceps Set

Reference	Description	Qty
MCEN54K*	Dessi Bipolar Forceps Set	1 Set
Set includes:		
MCEN54	Dessi Bipolar Forceps (6.0 mm Jaws) Working Length 115.0 mm (4.5")	1
MCEN540	Instrument Tray for Forceps and Cable 17.5 cm Wide x 28.5 cm Long x 5.0 cm Deep	1
CP393	Bifurcated (banana plug) Bipolar Cable (322.0 cm)	1

* Note: This set is only available in the US. If you are outside the US, please order the above items individually.



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

Manufacturer:



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France