

## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### Sataloff Laryngoscopes

Designed in conjunction with Robert T. Sataloff, MD

The Sataloff Laryngoscopes have a large proximal opening for binocular microscope vision. They offer dual light fibers and dual suction for maximum illumination and visualization during microsurgery. The elevated triangular tips have been blunted which may help reduce tissue trauma.


The anterior commissure laryngoscopes (MCL120, MCL121, MCL122, and MCL123) have been designed to help provide the ability to examine the vocal cords at the anterior commissure. In addition, the MCL125 anterior commissure laryngoscope was designed for patients with limited laryngeal exposure.

The Vallecula Laryngoscope (MCL124) has a wide proximal opening and a high blunted lift on the distal tip.




### Ordering Information


Product	Description	Quantity
MCL120	Laryngoscope with Dual Suction and Light Fibers, Small	1
MCL121	Laryngoscope with Dual Suction and Light Fibers, Medium	1
MCL122	Laryngoscope with Dual Suction and Light Fibers, Large	1
MCL123	Laryngoscope with Dual Suction and Light Fibers, Small-Medium	1
MCL124	Vallecula Laryngoscope with Dual Suction and Light Fibers	1
MCL125	Anterior Commissure Laryngoscope with Dual Suction and Light Fibers For patients with limited anterior exposure	1




MCL120




MCL123




MCL124



MCL121



MCL122



MCL125

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance S.A.S.

Le Pavillon, 03160

Saint-Aubin-Le Monial ▪ France