


# Integra<sup>®</sup> MicroFrance<sup>®</sup> Sataloff Instrument Set

Designed in conjunction with Robert T. Sataloff, MD, DMA, FACS

Delicate laryngeal microsurgery requires very sharp, precise, fine-tipped instruments to remove pathology without provoking scar formation. MicroFrance<sup>®</sup> Sataloff Instruments are designed to help the surgeon perform minimally traumatic ENT procedures.



# Sataloff Instrument Set



The unique heart-shaped graspers feature atraumatic jaws with high quality serrations, ideal for mini-microflap procedures.

## MCLS40\*

### Forceps



Downbiting Cup Forceps  
MCL-S18



Right Heart-Shaped Grasper  
MCLS21



Left Heart-Shaped Grasper  
MCLS22



Right Alligator Forceps  
MCL-S23



Left Alligator Forceps  
MCL-S24



Right Alligator Forceps, Fine  
MCL-S39



Left Alligator Forceps, Fine  
MCL-S41

### Injection Needles



Straight Injection Needle  
MCL-S4

- 27 gauge
- For use with MCL-S0 and MCL-S27
- Single-use, non-sterile, box of 6



Right Angle Injection Needle  
MCL-S5

- 27 gauge
- For use with MCL-S0 and MCL-S27
- Single-use, non-sterile, box of 6



Injection Needle Connector  
MCL-S0



Angled Adapter for MCL-S0  
MCL-S27



Percutaneous Injection Needle  
MCL-S35

- 18 gauge
- For use with MCL55

### Suctions



French Velvet Eye Suction  
MCL-S16

- 0.89mm diameter



Velvet Eye Suction with Retractor Plate Below  
MCL-S17

- 0.89mm diameter



Velvet Eye Suction with Retractor Plate Above  
MCLS-17-1

- 0.89mm diameter

### Miscellaneous



Cyst Curette  
MCL-S38



Instrument Tray  
MCLS28

- 25.5cm wide x 54.5cm long x 7.5cm deep
- Includes MCL-S36 and MCL-S37

### Universal Handle and Instruments

Universal Handle  
MCL-S19

- 2 per set
- For use with Universal Handle Instruments



Curved Spatula  
MCL-S1



Double Hook  
MCL-S2



Cotton Carrier  
MCL-S3



Straight Spatula  
MCL-S7



Fine Angle Spatula  
MCL-S8



Sharp Right Angle Hook  
MCL-S9



Ball Dissector  
MCL-S10



Oblique Blunt Ball Dissector  
MCL-S11



Short Right Angle Ball Dissector  
MCL-S12



Long Right Angle Ball Dissector  
MCL-S13



Small Mirror  
MCL-S14



Large Mirror  
MCL-S15

### Knife Handle and Tips

Knife Handle  
MCL-S30

- 2 per set
- For use with single-use Knife Tips



Vascular Knife  
MCL-S31

- Single-use, sterile, box of 6



Mini Micro Flap Knife  
MCL-S34

- Single-use, sterile, box of 6



Sharp Knife  
MCL-S32

- Single-use, sterile, box of 6



Sickle Knife  
MCL-S33

- Single-use, sterile, box of 6

\*This set is only available in the US. If you are outside the US, items in the set must be ordered individually.



### Sataloff Microlaryngeal Instrument Accessories (not included in the set)



#### Sharp Knife

##### MCL-S6

- 7mm blade, non-sterile, box of 6
- Overall length 252mm (9.9")
- For use with MCL-S19



#### Vascular Knife

##### MCL-S20

- Overall length 250mm (9.9")
- Non-sterile, 1 per box
- For use with MCL-S19



#### Sickle Knife

##### MCL-S25

- 7mm blade, non-sterile, box of 6
- Overall length 251mm (9.9")
- For use with MCL-S19



#### Curved Injection Needle

##### MCL-S42

- 22 gauge
- Single use, non-sterile, box of 6
- For use with MCL-S0 and MCL-S27



#### Straight Injection Needle

##### MCL-S43

- 22 gauge
- Single use, non-sterile, box of 6
- For use with MCL-S0 and MCL-S27

#### Instrument Tray with Lid

##### MCL-S36

- One level
- 22.5cm wide x 49.5cm long

#### Instrument Tray for Sataloff Picks, Handles and Injection System

##### MCL-S37

- One level
- 8.9 wide x 19.5 long
- Holds 16 picks, 2 handles, 12 needles

The MicroFrance brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)

Manufacturer:



Integra MicroFrance S.A.S.

Le Pavillon

03160 Saint-Aubin-Le Monial ▪ France